



Miltenyi Biotec

# CliniMACS® Technology in stem cell transplantation



CliniMACS Prodigy®

# Enrichment of CD34<sup>+</sup> stem cells

## Overcome the challenges of HSCT

The scope of hematopoietic stem cell transplantation (HSCT) keeps expanding from malignant to non-malignant diseases. However, severe side effects like graft-versus-host disease (GVHD), infections, and poor graft function remain challenges of the treatment. We developed our CliniMACS Systems for GMP-compliant cell depletion and enrichment strategies to help clinicians overcome these challenges.

### Passive T and B cell depletion

The immunomagnetic *in vitro* enrichment of CD34<sup>+</sup> cells, using the CliniMACS System, is a very potent and robust technology for T cell depletion. Studies have shown 10<sup>4</sup> to 10<sup>5</sup>-fold depletion rates<sup>1</sup>, leading to an effective GVHD prevention.<sup>2-5</sup> In the US, the CliniMACS CD34 Reagent System was registered as sole GVHD prophylaxis in allogeneic HSCT, from an HLA-identical sibling donor in adult patients with acute myeloid leukemia (AML) in first complete remission. This technique has been widely used for more than two decades and is continuously improving in terms of full automation on the CliniMACS Prodigy.<sup>6</sup> However, the use of CD34<sup>+</sup> cell enrichment as GVHD prevention strategy, may contribute to delayed immune reconstitution and higher non-relapse mortality (NRM).<sup>7</sup> Yet, the addition of memory T cells to the graft may solve this challenge. Memory T cells were shown to support the immune reconstitution and significantly improve the NRM rate, while maintaining excellent GVHD prevention.<sup>8</sup>



**Figure 1:** Connecting CliniMACS CD34 Reagents to the CliniMACS Prodigy.

### Stem cell boost (SCB) – rescue therapy for poor graft function (PGF)

PGF is a complication after allogeneic HSCT, occurring in 5–27% of the cases.<sup>9</sup> It is associated with increased morbidity and mortality due to severe infections, hemorrhagic complications, and organ failure caused by iron-overload. A boost of selected CD34<sup>+</sup> stem cells from the initial donor, fresh or cryopreserved, without any further conditioning and GVHD prophylaxis, has been described as an option to overcome persistent PGF.<sup>10</sup>

#### Poor graft function

Measures to diagnose PGF: two or three cytopenias >2 weeks, after day +28 in the presence of donor chimerism >5%

EBMT Handbook 2019<sup>20</sup>

As published, application of stem cell boost was applied independent of donor type and transplantation approach, also after failed previous therapy with growth factors. Clinical results show a 72–100% rapid and sustained lineage recovery with low rates of GVHD, even in high risk and mixed chimerism patients. Three years overall survival rates of 40–63% were reported.<sup>10-16</sup> In addition, immune reconstitution of all major lymphocyte populations was shown to be improved in 80% of the cases within four weeks.<sup>17</sup>

Other approaches like SCB as prevention of PGF after post-transplantation cyclophosphamide (PTCy) or treatment of mixed chimerism of secondary immunodeficiency (SID) patients after reduced intensity conditioning (RIC) HSCT are currently being investigated.<sup>18,19</sup>

Overall response	Complete response	aGVHD	cGVHD	Survival rate, median follow up 42 months
80%	72%	17%	18%	54%

**Table 1:** Efficacy and safety of CD34<sup>+</sup> selected SCB for PGF after allo-HSCT, from seven studies (N = 209) adapted from Shahzad, M. *et al.* (2021).<sup>27</sup>

## The haplo-cord approach

CD34<sup>+</sup> cell enrichment has also been applied in combined haplo-cord transplantations, in which a CD34<sup>+</sup> cell-enriched haploidentical transplant bridges the gap until the primary engraftment of the cord blood transplant takes place.<sup>21-23</sup>

As described, the idea of a combined haplo-cord transplantation is based on the observation that engraftment after haploidentical transplantation occurs faster than after cord blood transplantation. To provide the benefit of faster engraftment, the patient receives a bridging haploidentical transplant together with the cord blood at the same time. Clinical data show that with this approach even cord blood units containing very low stem cell numbers can be transplanted successfully and lead to sustained engraftment.<sup>21-23</sup>

## CD34<sup>+</sup> cell enrichment in autologous HSCT

CD34<sup>+</sup> cell enrichment was originally developed for passive tumor cell depletion from autologous stem cell grafts. Currently, the technique is used for some tumors of early childhood<sup>24,25</sup> or non-Hodgkin lymphoma and other diseases in adults.<sup>2,3</sup> Furthermore, CD34<sup>+</sup> cell enrichment has been used to deplete autoreactive cells from stem cell grafts for severe refractory autoimmune diseases, like systemic lupus erythematosus and systemic sclerosis.<sup>4,5,26</sup>



# Active T and B cell depletion

## Prevention of GVHD


After active T and B cell depletion, the graft contains CD34<sup>+</sup> stem cells, CD34<sup>-</sup> stem cells, other progenitor cells, natural killer (NK) cells, and other members of the innate immune system that might have engraftment facilitating effects. Several strategies for depletion of distinct T cell subsets have been developed, such as the CliniMACS TCR α/β/CD19 Depletion System.<sup>28,29</sup>

### Active depletion – TCRα/β/CD19 depletion

The process results in the depletion of alloreactive TCRα/β<sup>+</sup> T cells. At this, various cell populations like stem cells and immune effectors cells are retained in the cellular product. Immune effector cells, such as NK cells and the TCRγ/δ<sup>+</sup> T cells, are reported to induce GVL/T effects while the potential risk of inducing GVHD may be reduced.<sup>30,31</sup>

### TCRγ/δ<sup>+</sup> T cells – a fascinating cell type

Their unique set of functions include the ability to directly lyse infected or stressed cells, the production of cytokines and chemokines, and antigen presentation comparable to dendritic cells. Furthermore, they are thought to be highly effective against tumor cells and common infections.



**Figure 2: T cell receptors of TCRα/β<sup>+</sup> and TCRγ/δ<sup>+</sup> T cells.**

One major difference lies in the T cell receptor which is composed of different chains, α and β or γ and δ. TCRγ/δ<sup>+</sup> T cells are not activated by MHC-presented antigens, which might contribute to the low alloreactivity of TCRγ/δ<sup>+</sup> T cells. These cells become activated by structures of bacterial walls or heat shock proteins, for example.

### Results from clinical applications

First clinical results using TCRα/β or TCRα/β/CD19 depleted stem cell grafts were published in 2011 by Lang *et al.*, from the Children’s Hospital University Tübingen, Germany.<sup>32</sup> Since then, the number of publications has increased tremendously. The technique has been applied in:

- malignant and non-malignant diseases
- myeloablative and reduced intensity conditioning regimens
- pediatric and adult patients
- settings with reduced or no post-transplant immune suppression
- haploidentical and matched unrelated HSCT

(see references 30–40)

	Rome <sup>33</sup>	Newcastle <sup>34</sup>	Utrecht <sup>35</sup>
CD34 <sup>+</sup> stem cells ×10 <sup>6</sup> /kg	15.8	17.8	6.1
TCRαβ <sup>+</sup> T cells ×10 <sup>3</sup> /kg	40	3.3	20
TCRγδ <sup>+</sup> T cells ×10 <sup>6</sup> /kg	9.4	n/a	5.1
NK cells ×10 <sup>6</sup> /kg	38.2	5.15	18.4
CD20 <sup>+</sup> B cells ×10 <sup>4</sup> /kg	4	4	41.3
No. of patients	23	25	35

**Table 2:** Graft composition of haploidentical HSCT from three different sites.

The observed rates for primary engraftment are described as ‘fast’ in various publications, with ten to 16 days for neutrophil recovery.<sup>33,36–40</sup> Also, the immune reconstitution data are rated as ‘remarkably fast’ by Balashov *et al.*<sup>38</sup> Bertaina *et al.*<sup>33</sup> published the clinical results from 23 children with non-malignant disorders, receiving haploidentical HSCT and reported a two-year probability of disease-free survival of 91.1%. In a study with an adult patient population suffering from high-risk leukemias, it was shown that infusion of grafts depleted of TCRα/β<sup>+</sup>/CD19<sup>+</sup> cells was safe and effective, resulting in rapid donor hematopoietic engraftment and early expansion of donor-derived T lymphocytes.<sup>37</sup>

High rates of disease-free survival, overall survival, and GVHD-free-relapse-free survival rates were observed in the company initiated multicenter prospective trial. Included were heavily pretreated adult and pediatric patients with malignant and non-malignant diseases using reduced-intensity conditioning regimen. The results showed a fast reconstitution of TCR $\gamma\delta^+$  T cells and NK cells in the early posttransplant period, which translated into a favorable incidence of infectious complications. The low rate of 10% acute graft-versus-host disease (aGVHD) and no cases of grade III–IV aGVHD impressively highlights the strength of *ex vivo* TCR $\alpha/\beta^+$ /CD19 $^+$  depletion.<sup>40</sup>



# Manufacturing of memory T cell products

## A versatile concept for transferring immunity

Memory T cells provide immunity against many viral infections that the donors have experienced during their lifetime. In contrast to naive T cells, memory T cells are less alloreactive in an allogeneic HSCT setting.<sup>41-43</sup> CliniMACS CD45RA-depleted cell products have been used to treat and prevent opportunistic infections and to enhance immune reconstitution in the setting of allogeneic HSCT.<sup>43-45</sup>

CD45RA is an isoform of the CD45 family. It is present on naive T cells, on other leukocytes, and on some parts of hematopoietic stem cells in mobilized leukapheresis products. CD45RA is a proven marker to deplete alloreactivity in one step while central and effector memory cells are preserved.<sup>41,43,46-48</sup>

Log depletion	Leukapheresis product
4.7 (naive T cells)	(Bleakley, 2014)
4.4 (3.4-4.7; CD45RA <sup>+</sup> cells)	(Teschner, 2013)
3.6 (2.3-3.95; CD45RA <sup>+</sup> cells)	(Triplett, 2015)

**Table 3:** Log depletions obtained by CliniMACS Plus CD45RA Depletion of mobilized leukapheresis products.<sup>41,42,49</sup>

### Graft engineering approach on an *ex vivo* T cell depletion platform

The intention of an *ex vivo* T cell depletion strategy is to provide grafts consisting of stem cells, combined with memory T cells, depleted of alloreactive naive T cell. The CliniMACS Platform allows for various combinations.

- The combination of CD34<sup>+</sup> cell enrichment and subsequent CD45RA<sup>+</sup> cell depletion from the CD34-negative cell fraction of the same mobilized leukapheresis products, has first been described by Marie Bleakley *et al.*<sup>41</sup> This approach has been used in different transplantation settings, as matched related and haploidentical transplantations (Seattle, Madrid, Paris).<sup>8,43,50,51</sup>
- Researcher from St. Judes Children's Hospital combined a CD34 enriched graft with CD45RA<sup>+</sup> depleted memory T cells that were obtained from an additional leukapheresis product. In this treatment protocol, CD56<sup>+</sup> enriched NK cells were added as well.<sup>52,53</sup>
- Other researchers developed strategies in which TCR $\alpha/\beta$ -depleted grafts were combined with CD45RA-depleted memory T cells. This approach was mainly used in haploidentical transplantation settings.<sup>54</sup>

## Donor Lymphocyte Infusion (DLI) – memory T cells posttransplant

DLIs, depleted of CD45RA<sup>+</sup> cells, have been infused posttransplant as infection prophylaxis. This approach might be applied in combination with any kind of GVHD prophylaxis, such as *ex vivo* or *in vivo* T cell depletion strategies.

- The children's hospital in Moscow, Russia, reported on prophylactic infusions of memory T cell DLI, depleted of CD45RA<sup>+</sup> cells, after receiving a TCR $\alpha$ / $\beta$ -depleted haploidentical transplant. The patients received up to three infusions in monthly intervals.<sup>44</sup>
- Prophylactic infusions of memory T cells for patients receiving T cell-depleted grafts from identical sibling donors have also been reported for the Duke University Medical Center, US.<sup>55</sup>
- In a study from Milan, Italy, patients received prophylactic memory T cell infusions after treatment with a haploidentical transplant that was based on a post-cyclophosphamide based GVHD prophylaxis.<sup>56</sup>

The use of memory T cells to treat active drug-resistant infections was reported too. The treatment of active infections after haploidentical HSCT in steroid-refractory cases has been demonstrated in several case reports. In these cases, the memory T cell products were obtained from original stem cell donors.<sup>45,57</sup>

## Memory T cells in COVID-19


It was demonstrated that donor cells in a non-transplant setting can be helpful to overcome infections. Antonio Perez and his collaborators developed a method for delivering cell therapy products to COVID-19 patients. Using CliniMACS CD45RA Depletion, cell products were generated that are stored within a COVID-19 biobank from convalescent donors.<sup>58</sup> The clinical use of these cells has been reported within a phase I clinical trial.<sup>59</sup>



# Enrichment of antigen-specific T cells

## Enable anti-viral immunity

Antigen-specific T cells can be enriched using the CliniMACS Cytokine Capture System (CCS) (IFN-gamma) Technology. The CliniMACS Prodigy CCS (IFN-gamma) System allows for a fully automated specific labeling and enrichment procedure for antigen-specific IFN- $\gamma$  secreting CD4<sup>+</sup> and CD8<sup>+</sup> T cells. The specificity of enriched T cells can be determined by choosing the respective MACS<sup>®</sup> GMP PepTivator<sup>®</sup> Peptide Pools. A growing repertoire of PepTivators for various viral antigens like AdV, BKV, EBV, and HCMV is available. Additionally, by using MACS GMP PepTivator NY-ESO-1 or WT1, cancer antigen-specific T cells may also be isolated.



**Figure 3:** Overview of the automated workflow steps of the CliniMACS Prodigy CCS (IFN-gamma) System. Virus-specific T cells can be enriched based on their secretion of IFN- $\gamma$  after restimulation with appropriate antigens.

## Treatment of drug-resistant infections after allogeneic HSCT

Published clinical data are available for patients suffering from CMV, AdV, BKV, and EBV infections after receiving allogeneic HSCT. In these cases the infections could not be resolved by standard antiviral treatment. Initially, cells were obtained from the original stem cell donor.<sup>60-63</sup> Today, an increasing number of reports show the possibility of using third-party donor-derived cell products that can provide different types of biobanks.<sup>64,65</sup> Currently a multi-national clinical phase III trial is investigating efficacy and safety of adoptive T cell transfer in immunocompromised individuals.<sup>66</sup>

## Treatment of patients after solid organ transplantation

Based on the experiences with HSCT patients, new strategies have been developed, using third party donor-derived cell products for patients post solid organ transplantation. Case reports are available for EBV-related lymphoma that has been treated with CCS derived antigen-specific T cells.<sup>67,68</sup>

## SARS-CoV-2 specific T cells against COVID-19 disease

Several groups have developed manufacturing processes for cell products from convalescent donors, targeting SARS-CoV-2 infections. Wing Leung and his co-workers describe a process for flexible and rapid manufacturing of virus-specific T cells (VST) once a donor is available.<sup>69</sup> Cooper *et al.* describes a process that starts with the enrichment of VSTs using the CliniMACS Cytokine Capture System. A culturing step is added to increase the harvest of VST. The final products can be provided within a biobank.<sup>70</sup> First phase I clinical trials have been started.<sup>71,72</sup>





 PepTivator® CEF M  
Class I Plus - prem  
grade 6 nmol/pept  
FOR RESEARCH USE ON  
Order no.: 130-098-425  
Contains stabilizer  
Store at 20 °C

PepTivator Peptide Pools are optimal for effective stimulation of antigen-specific T cells.

# CliniMACS® Plus System

## Automated cell separation for clinical-scale cell enrichment or depletion

The CliniMACS Plus System automates cell separation for clinical-scale enrichment of target cells or depletion of unwanted cells from blood products. Cell separation occurs in a closed and sterile system. A single-use tubing set with its integrated separation column enables the instrument to separate target cells from unwanted cells and collect the fractions in different bags.

### Sterile connections

Cellular starting material and buffers are connected for cell processing.

### Ease of use

The software guides the user through interactions.

### Flow of liquids

A peristaltic pump takes care for the precise flow of liquids.

### Automated cell separation

The magnet unit houses the separation column for immunomagnetic cell separation

### Controlled pathways

Pinch valves automatically close and open for a controlled fluid pathway within the tubing sets.




## Key features of the CliniMACS Plus System

- Clinical-scale cell separation
- Automated separation procedure
- Compatible for use in a GMP setting
- Closed system, using CE-marked medical devices
- Enrichment or depletion of dedicated cell populations
- Reproducible high purities, excellent yields and profound depletion efficiencies
- Established platform for *ex vivo* T cell depletion since 1997

## Modular platform

- Variable cell sources
- Wide range of applications
- Versatile separation strategies



**Figure 5:** Workflow on the CliniMACS Plus System. Target cells are always collected in the left-most bag.

# CliniMACS Prodigy® Platform

## Master the complexity of cell processing

The CliniMACS Prodigy integrates all cell processing steps, including sample preparation, cell washing, density gradient centrifugation, magnetic cell separation, cell activation, genetic modification, cell culture, and cell product formulation. The fully automated, sensor-controlled processes provide a high level of standardization and reproducibility. Hands-on time is reduced substantially.

As all steps are performed in single-use, closed and sterile tubing sets, the instrument also reduces cleanroom requirements significantly. In combination with the wide variety of CE-certified medical devices and MACS GMP Products manufactured by Miltenyi Biotec, the CliniMACS Prodigy facilitates the implementation of GMP-compliant cell processing.

### Network connectivity

The instrument is equipped to be connected to your local area network for process surveillance and data export.

### Flow of liquids


The peristaltic pump directs accurate volumes of liquids through the tubing set.

### Centrifugation and cell processing

The CentriCult™ Unit is a multi-purpose entity for cell processing and cell cultivation.

### Controlled pathways

Pinch valves automatically close and open for a controlled fluid pathway within the tubing set.





### Safe and sterile connections

Reagents, buffers, media and other materials can be easily connected via sterile spike ports or sterile welding.

### Easy to operate

The user-friendly software guides the operator through interactions. Predefined processes and the flexible programming suite allow for convenient and tailored cell manufacturing.

### Automated cell selection

The magnet unit houses the separation column for immunomagnetic cell separation.

### Effortless documentation

CliniMACS Materials are scanned for identification and documentation with the bar code reader.

### Sterile sealing

The MACS Tube Sealer is used for sealing and sterile removal of tubing components such as QC pouches.

# Selected references

1. Panch, S. R. *et al.* (2019) Robust selections of various hematopoietic cell fractions on the CliniMACS Plus Instrument. *Clin. Hematol. Int.* 1(3): 161–167.  
<https://pubmed.ncbi.nlm.nih.gov/34595426/>
2. Yahng, S. A. *et al.* (2014) Influence of *ex vivo* purging with CliniMACS CD34(+) selection on outcome after autologous stem cell transplantation in non-Hodgkin lymphoma. *Br. J. Haematol.* 164: 555–564.  
<https://pubmed.ncbi.nlm.nih.gov/24266323/>
3. Berger, M. D. *et al.* (2015) CD34+ selected versus unselected autologous stem cell transplantation in patients with advanced-stage mantle cell and diffuse large B-cell lymphoma. *Leuk. Res.* 39: 561–567.  
<https://doi.org/10.1016/j.leukres.2015.03.004>
4. Alchi, B. *et al.* (2013) Autologous haematopoietic stem cell transplantation for systemic lupus erythematosus: data from the European Group for Blood and Marrow Transplantation registry. *Lupus* 22: 245–253.  
<https://doi.org/10.1177/0961203312470729>
5. Henes, J. C. *et al.* (2012) Optimization of autologous stem cell transplantation for systemic sclerosis — a single-center longterm experience in 26 Patients with severe organ manifestations. *J. Rheumatol.* 39: 269–275.  
<https://doi.org/10.3899/jrheum.110868>
6. Avelilla, S. T. *et al.* (2016) How do I perform hematopoietic progenitor cell selection? *Transfusion* 56(5): 1008–1012.  
<https://doi.org/10.1111/trf.13534>
7. Luznik, L. *et al.* (2021) Randomized phase III BMT CTN trial of calcineurin inhibitor-free chronic graft-versus-host disease interventions in myeloablative hematopoietic cell transplantation for hematologic malignancies. *J. Clin. Oncol.* JCO2102293.  
<https://doi.org/10.1200/jco.21.02293>
8. Bleakley, M. *et al.* (2022) Naive t-cell depletion to prevent chronic graft-versus-host disease. *J. Clin. Oncol.* JCO2101755.  
<https://doi.org/10.1200/JCO.21.01755>
9. Lee, K. H. *et al.* (2004) Failure of trilineage blood cell reconstitution after initial neutrophil engraftment in patients undergoing allogeneic hematopoietic cell transplantation - Frequency and outcomes. *Bone Marrow Transplant.* 33(7): 729–734.  
<https://doi.org/10.1038/sj.bmt.1704428>
10. Ghobadi, A. *et al.* (2017) Fresh or cryopreserved CD34+-selected mobilized peripheral blood stem and progenitor cells for the treatment of poor graft function after allogeneic hematopoietic cell transplantation. *Biol. Blood Marrow Transplant.* 23(7): 1072–1077.  
<https://doi.org/10.1016/j.bbmt.2017.03.019>
11. Stasia, A. *et al.* (2014) CD34 selected cells for the treatment of poor graft function after allogeneic stem cell transplantation. *Biol. Blood Marrow Transplant.* 20: 1440–1443.  
<https://doi.org/10.1016/j.bbmt.2014.05.016>
12. Klyuchnikov, E. *et al.* (2014) CD34+-selected stem cell boost without further conditioning for poor graft function after allogeneic stem cell transplantation in patients with hematological malignancies. *Biol. Blood Marrow Transplant.* 20: 382–386.  
<https://doi.org/10.1016/j.bbmt.2013.11.034>
13. Haen, S.P. *et al.* (2015) Poor graft function can be durably and safely improved by CD34+-selected stem cell boosts after allogeneic unrelated matched or mismatched hematopoietic cell transplantation. *J. Cancer Res. Clin. Oncol.* 141: 2241–2251.  
<https://doi.org/10.1007/s00432-015-2027-x>
14. Askaa, B. *et al.* (2014) Treatment of poor graft function after allogeneic hematopoietic cell transplantation with a booster of CD34-selected cells infused without conditioning. *Bone Marrow Transplant.* 49: 720–721.  
<https://doi.org/10.1038/bmt.2014.5>
15. Chandra, S. *et al.* (2018) Post-transplant CD34+ selected stem cell “boost” for mixed chimerism after reduced-intensity conditioning hematopoietic stem cell transplantation in children and young adults with primary immune deficiencies. *Biol. Blood Marrow Transplant.* 24: 1527–1529.  
<https://doi.org/10.1016/j.bbmt.2018.03.013>
16. Mohy, R. *et al.* (2019) CD34+-selected stem cell “boost” for poor graft function after allogeneic hematopoietic stem cell transplantation. *Curr. Res. Transl. Med.* 67: 112–114.  
<https://doi.org/10.1016/j.retram.2018.12.003>
17. Cuadrado, M. *et al.* (2020) Predictors of recovery following allogeneic CD34+-selected cell infusion without conditioning to correct poor graft function. *Haematologica* 105: 2639–2646.  
<https://doi.org/10.3324/haematol.2019.226340>
18. Mainardi, C. *et al.* (2018) CD34+ selected stem cell boosts can improve poor graft function after paediatric allogeneic stem cell transplantation. *Br. J. Haematol.* 180: 90–99.  
<https://doi.org/10.1111/bjh.15012>
19. Abboud, R. *et al.* (2021) Can planned CD34+ stem cell boost prevent poor graft function after peripheral blood haploidentical hematopoietic transplantation? *Leuk. Lymphoma* 62: 749–751.  
<https://doi.org/10.1080/10428194.2020.1839657>
20. Carreras, E. *et al.* (2019) EBMT Handbook 2019. Hematopoietic stem cell transplantation and cellular therapies. 311.  
<https://doi.org/10.1007/978-3-030-02278-5>
21. Liu, H. & van Besien, K. (2015) Alternative donor transplantation – “mixing and matching”: the role of combined cord blood and haplo-identical donor transplantation (haplo-cord SCT) as a treatment strategy for patients lacking standard donors? *Curr. Hematol. Malig. Rep.* 10: 1–7.  
<https://doi.org/10.1007/s11899-014-0245-y>
22. Kwon, M. *et al.* (2014) Haplo-cord transplantation using CD34+ cells from a third-party donor to speed engraftment in high-risk patients with hematologic disorders. *Biol. Blood Marrow Transplant.* 20: 2015–2022.  
<https://doi.org/10.1016/j.bbmt.2014.08.024>
23. Van Besien, K. *et al.* (2019) Haploidentical vs haplo-cord transplant in adults under 60 years receiving fludarabine and melphalan conditioning. *Blood Adv.* 3: 1858–1867.  
<https://doi.org/10.1182/bloodadvances.2019000200>
24. Ballova, V. *et al.* (2008) Autologous stem cell transplantation with selected CD34+ cells and unmanipulated peripheral blood stem cells in patients with relapsed and refractory Hodgkin's lymphoma: a single centre experience. *Neoplasma* 55: 428–436.  
<https://pubmed.ncbi.nlm.nih.gov/18665754/>
25. Marabelle, A. *et al.* (2011) CD34+ immunoselection of autologous grafts for the treatment of high-risk neuroblastoma. *Pediatr. Blood Cancer* 56: 134–142.  
<https://doi.org/10.1002/pbc.22840>
26. van Laar, J. M. *et al.* (2014) Autologous hematopoietic stem cell transplantation vs intravenous pulse cyclophosphamide in diffuse cutaneous systemic sclerosis. A randomized clinical trial. *JAMA* 311: 2490–2498.  
<https://doi.org/10.1001/jama.2014.6368>
27. Shahzad, M. *et al.* (2021) Outcomes with CD34-selected stem cell boost for poor graft function after allogeneic hematopoietic stem cell transplantation: a systematic review and meta-analysis. *Transplant. Cell. Ther.* 27: 877.e1–877.e8.  
<https://doi.org/10.1016/j.jctc.2021.07.012>
28. Chaleff, S. *et al.* (2007) A large-scale method for the selective depletion of  $\alpha\beta$  T lymphocytes from PBSC for allogeneic transplantation. *Cytotherapy* 9:746–754.  
<https://doi.org/10.1080/14653240701644000>
29. Schumm, M. *et al.* (2013) Depletion of T-cell receptor alpha/beta and CD19 positive cells from apheresis products with the CliniMACS device. *Cytotherapy* 15: 1253–1258.  
<https://doi.org/10.1016/j.jcyt.2013.05.014>
30. Locatelli, F. *et al.* (2013) At the bedside: innate immunity as an immunotherapy tool for hematological malignancies. *J. Leukoc. Biol.* 94: 1141–1157.  
<https://doi.org/10.1189/jlb.0613343>
31. Vantourout, P. *et al.* (2013) Six-of-the-best: unique contributions of  $\gamma\delta$  T cells to immunology. *Nat. Rev. Immunol.* 13: 88–100.  
<https://doi.org/10.1038/nri3384>
32. Lang, P. *et al.* (2011) First clinical results with  $\alpha\beta$ + T-cell depleted haploidentical stem cells in children. *Bone Marrow Transplant.* 46: P558–S135.  
<https://doi.org/10.1038/bmt.2011.47>
33. Bertaina, A. *et al.* (2014) HLA-haploidentical stem cell transplantation after removal of  $\alpha\beta$ + T and B cells in children with nonmalignant disorders. *Blood* 124: 822–826.  
<https://doi.org/10.1182/blood-2014-03-563817>
34. Shah, R. *et al.* (2017) T-cell receptor  $\alpha\beta$ + and CD19+ cell-depleted haploidentical and mismatched hematopoietic stem cell transplantation in primary immune deficiency. *J. Allergy Clin. Immunol.* 141: 1417–1426.e1.  
<https://doi.org/10.1016/j.jaci.2017.07.008>
35. De Witte, M. A. *et al.* (2021)  $\alpha\beta$  T-cell graft depletion for allogeneic HSCT in adults with hematological malignancies. *Blood Adv.* 5: 240–249.  
<https://doi.org/10.1182/bloodadvances.2020002444>
36. Lang, P. *et al.* (2015) Improved immune recovery after transplantation of TCR $\alpha\beta$ /CD19-depleted allografts from haploidentical donors in pediatric patients. *Bone Marrow Transplant.* 50: 6–10.  
<https://doi.org/10.1038/bmt.2015.87>
37. Prezioso, L. *et al.* (2019) Haploidentical hematopoietic stem cell transplantation in adults using the  $\alpha\beta$ TCR/CD19-based depletion of G-CSF-mobilized peripheral blood progenitor cells. *Bone Marrow Transplant.* 54: 698–702.  
<https://doi.org/10.1038/s41409-019-0608-z>

38. Balashov, D. *et al.* (2015) Single-center experience of unrelated and haploidentical stem cell transplantation with TCR $\alpha\beta$  and CD19 depletion in children with primary immunodeficiency syndromes. *Biol. Blood Marrow Transplant.* 21: 1955–1962. <https://doi.org/10.1016/j.bbmt.2015.07.008>
39. Maschan, M. *et al.* (2016) TCR-alpha/beta and CD19 depletion and treosulfan-based conditioning regimen in unrelated and haploidentical transplantation in children with acute myeloid leukemia. *Bone Marrow Transplant.* 51: 668–674. <https://doi.org/10.1038/bmt.2015.343>
40. Bethge, W. A. *et al.* (2021) Results of a multicenter phase I/II trial of TCR $\alpha\beta$  and CD19-depleted haploidentical hematopoietic stem cell transplantation for adult and pediatric patients. *Bone Marrow Transplant.* Online ahead of print <https://doi.org/10.1038/s41409-021-01551-z>
41. Bleakley, M. *et al.* (2014) Engineering human peripheral blood stem cell grafts that are depleted of naive T cells and retain functional pathogen-specific memory T cells. *Biol. Blood Marrow Transplant.* 20: 705–716. <https://doi.org/10.1016/j.bbmt.2014.01.032>
42. Teschner, D. *et al.* (2014) Depletion of naive T cells using clinical grade magnetic CD45RA beads: a new approach for GVHD prophylaxis. *Bone Marrow Transplant.* 49: 138–144. <https://doi.org/10.1038/bmt.2013.114>
43. Bleakley, M. *et al.* (2015) Outcomes of acute leukemia patients transplanted with naive T cell–depleted stem cell grafts. *J. Clin. Invest.* 125: 2677–2689. <https://doi.org/10.1172/JCI81229>
44. Dunaikina, M. *et al.* (2021) Safety and efficacy of the low-dose memory (CD45RA-depleted) donor lymphocyte infusion in recipients of  $\alpha\beta$  T cell-depleted haploidentical grafts: results of a prospective randomized trial in high-risk childhood leukemia. *Bone Marrow Transplant.* 56: 1614–1624. <https://doi.org/10.1038/s41409-021-01232-x>
45. Brodzski, N. *et al.* (2016) Novel treatment of severe combined immunodeficiency utilizing *ex vivo* T-cell depleted haploidentical hematopoietic stem cell transplantation and CD45RA<sup>+</sup> depleted donor lymphocyte infusions. *Orphanet J. Rare Dis.* 11: 5. <https://doi.org/10.1186/s13023-016-0385-3>
46. Fritsch, G. *et al.* (1993) Rapid discrimination of early CD34<sup>+</sup> myeloid progenitors using CD45-RA analysis. *Blood* 81: 2301–2309. <https://doi.org/10.1182/blood.V81.9.2301.2301>
47. Distler, E. *et al.* (2011) Alloreactive and leukemia-reactive T cells are preferentially derived from naive precursors in healthy donors: implications for immunotherapy with memory T cells. *Haematologica* 96: 1024–1032. <https://doi.org/10.3324%2Fhaematol.2010.037481>
48. Bremm, M. *et al.* (2021) Depletion of CD45RA<sup>+</sup> T cells: Advantages and disadvantages of different purification methods. *Journal of Immunological Methods* 492: 112960. <https://doi.org/10.1016/j.jim.2021.112960>
49. Triplett, B. M. *et al.* (2015) Rapid memory T-cell reconstitution recapitulating CD45RA-depleted haploidentical transplant graft content in patients with hematologic malignancies. *Bone Marrow Transplant* 50: 968–977. <https://doi.org/10.1038/bmt.2014.324>
50. Touzot, F. *et al.* (2015) CD45RA depletion in HLA-mismatched allogeneic hematopoietic stem cell transplantation for primary combined immunodeficiency: a preliminary study. *J. Allergy Clin. Immunol.* 135: 1303–1309. <https://doi.org/10.1016/j.jaci.2014.08.019>
51. Sisinni, L. *et al.* (2018) Unexpected high incidence of human herpesvirus-6 encephalitis after naive T cell–depleted graft of haploidentical stem cell transplantation in pediatric patients. *Biol. Blood Marrow Transplant.* 24: 2316–2323. <https://doi.org/10.1016/j.bbmt.2018.07.016>
52. Triplett, B. M. *et al.* (2018) Selective T-cell depletion targeting CD45RA reduces viremia and enhances early T-cell recovery compared with CD3-targeted T-cell depletion. *Transpl. Infect. Dis.* 2018: e12823. <https://doi.org/10.1111/tid.12823>
53. Mamcarz, E. *et al.* (2020) Improved survival rate in T-cell depleted haploidentical hematopoietic cell transplantation over the last 15 years at a single institution. *Bone Marrow Transplant* 55: 929–938. <https://doi.org/10.1038/s41409-019-0750-7>
54. Wu, I. Q. *et al.* (2021) HLA-Haploidentical hematopoietic cell transplantation of *ex vivo* TCR $\alpha\beta$ -depleted grafts with CD45RA-depleted memory T cell add-back in adults and children with hematological malignancies: 4-year follow-up of multicenter study in Singapore. *Blood* 138 (suppl. 1): 555. <https://doi.org/10.1182/blood-2021-152372>
55. Maung, K. K. *et al.* (2021) Phase I dose escalation study of naive T-cell depleted donor lymphocyte infusion following allogeneic stem cell transplantation. *Bone Marrow Transplant* 56: 137–143. <https://doi.org/10.1038/s41409-020-0991-5>
56. Castagna, L. *et al.* (2021) Feasibility and efficacy of CD45RA<sup>+</sup> depleted donor lymphocytes infusion after haploidentical transplantation with post-transplantation cyclophosphamide in patients with hematological malignancies. *Transplantation and Cellular Therapy* 27: 478.e1–478.e5. <https://doi.org/10.1016/j.jtct.2021.03.010>
57. Park, H. J. *et al.* (2020) Successful treatment of refractory CMV colitis after haploidentical HSCT with post-transplant cyclophosphamide using CD45RA<sup>+</sup> depleted donor lymphocyte infusion. *Bone Marrow Transplant* 55: 1674–1676. <https://doi.org/10.1038/s41409-019-0685-z>
58. Ferreras, C. *et al.* (2021) SARS-CoV-2-specific memory T lymphocytes from COVID-19 convalescent donors: identification, biobanking, and large-scale production for adoptive cell therapy. *Front. Cell Dev. Biol.* 9: 620730. <https://doi.org/10.3389/fcell.2021.620730>
59. Perez-Martinez, A. *et al.* (2021) Phase I dose-escalation single centre clinical trial to evaluate the safety of infusion of memory T cells as adoptive therapy in COVID-19 (RELEASE). *EclinicalMedicine* 39: 101086. <https://doi.org/10.1016/j.eclinm.2021.101086>
60. Feuchtinger, T. *et al.* (2010) Adoptive transfer of pp65-specific T cells for the treatment of chemorefractory cytomegalovirus disease or reactivation after haploidentical and matched unrelated stem cell transplantation. *Blood* 116: 4360–4367. <https://doi.org/10.1182/blood-2010-01-262089>
61. Meiji, P. *et al.* (2012) Effective treatment of refractory CMV reactivation after allogeneic stem cell transplantation with in vitro-generated CMV pp65-specific CD8<sup>+</sup> T-cell lines. *J. Immunother.* 35: 621–628. <https://doi.org/10.1097/CJI.0b013e31826e35f6>
62. Icheva, V. *et al.* (2012) Adoptive transfer of Epstein-Barr virus (EBV) nuclear antigen 1-specific T cells as treatment for EBV reactivation and lymphoproliferative disorders after allogeneic stem-cell transplantation. *J. Clin. Oncol.* 31: 39–48. <https://doi.org/10.1200/JCO.2011.39.8495>
63. Pello, O. M. *et al.* (2017) BKV-specific T cells in the treatment of severe refractory haemorrhagic cystitis after HLA-haploidentical haematopoietic cell transplantation. *Eur. J. Haem.* 98: 632–634 <https://doi.org/10.1111/ejh.12848>
64. Qian, C. *et al.* (2017) Curative or pre-emptive adenovirus-specific T cell transfer from matched unrelated or third party haploidentical donors after HSCT, including UCB transplantations: a successful phase I/II multicenter clinical trial. *Journal of Hematology & Oncology* 10: 102. <https://doi.org/10.1186/s13045-017-0469-0>
65. Kállay, K. *et al.* (2018) Early experience with CliniMACS Prodigy CCS (IFN-gamma) system in selection of virus-specific T cells from third-party donors for pediatric patients with severe viral infections after hematopoietic stem cell transplantation. *J. Immunother.* 41: 158–163. <https://doi.org/10.1097/CJI.000000000000197>
66. Kaeuferle, T. *et al.* (2019) Strategies of adoptive T-cell transfer to treat refractory viral infections post allogeneic stem cell transplantation. *Journal of Hematology & Oncology* 12: 13. <https://doi.org/10.1186/s13045-019-0701-1>
67. Hansen, B. T. *et al.* (2021) Adoptive cell transfer of allogeneic Epstein-Barr virus-specific T lymphocytes for treatment of refractory EBV-associated posttransplant smooth muscle tumors: a case report. *Front. Immunol.* 12: 727814. <https://doi.org/10.3389/fimmu.2021.727814>
68. Schultze-Florey, R. E. *et al.* (2018) Dissecting Epstein-Barr virus-specific T-cell responses after allogeneic EBV-specific T-cell transfer for central nervous system posttransplant lymphoproliferative disease. *Front. Immunol.* 9: 1475. <https://doi.org/10.3389/fimmu.2018.01475>
69. Leung, W. *et al.* (2020) Rapid production of clinical-grade SARS-CoV-2-specific T cells. *Adv. Cell Gene. Ther.* 2020; 00: e101 <https://doi.org/10.1002/acg2.101>
70. Cooper, R. S. *et al.* (2021) Rapid GMP-compliant expansion of SARS-CoV-2-specific T cells from convalescent donors for use as an allogeneic cell therapy for COVID-19. *Front. Immunol.* 11: 598402 <https://doi.org/10.3389/fimmu.2020.598402>
71. clinicaltrials.gov: [NCT04324463](https://clinicaltrials.gov/ct2/show/study/NCT04324463)
72. clinicaltrials.gov: [NCT04742595](https://clinicaltrials.gov/ct2/show/study/NCT04742595)



**Miltényi Biotec**

**Germany/Austria**

Miltényi Biotec B.V. & Co. KG  
Friedrich-Ebert-Straße 68  
51429 Bergisch Gladbach  
Germany  
Phone +49 2204 8306-0  
Fax +49 2204 85197  
macsde@miltényi.com

**USA/Canada**

Miltényi Biotec, Inc.  
2303 Lindbergh Street  
Auburn, CA 95602, USA  
Phone 800 FOR MACS  
Phone +1 866 811 4466  
Fax +1 877 591 1060  
macsus@miltényi.com

**Australia**

Miltényi Biotec  
Australia Pty. Ltd.  
Unit 11, 2 Eden Park Drive  
Macquarie Park, NSW 2113  
Australia  
Phone +61 2 8877 7400  
Fax +61 2 9889 5044  
macsau@miltényi.com

**Benelux**

Miltényi Biotec B.V.  
Sandifortdreef 17  
2333 ZZ Leiden  
The Netherlands  
macsnl@miltényi.com

**Customer service**

**The Netherlands**  
Phone 0800 4020120  
Fax 0800 4020100

**Customer service Belgium**

Phone 0800 94016  
Fax 0800 99626

**Customer service Luxembourg**

Phone 800 24971  
Fax 800 24984

**China**

Miltényi Biotec Technology &  
Trading (Shanghai) Co., Ltd.  
Room 401  
No. 1077, Zhangheng Road  
Pudong New Area  
201203 Shanghai, P.R. China  
Phone +86 21 6235 1005-0  
Fax +86 21 6235 0953  
macscn@miltényi.com

**France**

Miltényi Biotec SAS  
10 rue Mercoeur  
75011 Paris, France  
Phone +33 1 56 98 16 16  
macsfr@miltényi.com

**Hong Kong**

Miltényi Biotec Hong Kong Ltd.  
Unit 301, Lakeside 1  
No. 8 Science Park West Avenue  
Hong Kong Science Park  
Pak Shek Kok, New Territories  
Hong Kong  
Phone +852 3751 6698  
Fax +852 3619 5772  
macshk@miltényi.com

**Italy**

Miltényi Biotec S.r.l.  
Via Paolo Nanni Costa, 30  
40133 Bologna  
Italy  
Phone +39 051 6 460 411  
Fax +39 051 6 460 499  
macsit@miltényi.com

**Japan**

Miltényi Biotec K.K.  
NEX-Eitai Building 5F  
16-10 Fuyuki, Koto-ku  
Tokyo 135-0041, Japan  
Phone +81 3 5646 8910  
Fax +81 3 5646 8911  
macsjp@miltényi.com

**Nordics and Baltics**

Miltényi Biotec Norden AB  
Medicon Village  
Scheeleorget 1  
223 81 Lund  
Sweden  
macsse@miltényi.com

**Customer service Sweden**

Phone 0200 111 800  
Fax +46 280 72 99

**Customer service Denmark**

Phone 80 20 30 10  
Fax +46 46 280 72 99

**Customer service**

**Norway, Finland, Iceland,  
and Baltic countries**  
Phone +46 46 280 72 80  
Fax +46 46 280 72 99

**Singapore**

Miltényi Biotec Asia Pacific Pte Ltd.  
438B Alexandra Road, Block B  
Alexandra Technopark  
#06-01  
Singapore 119968  
Phone +65 6238 8183  
Fax +65 6238 0302  
macssg@miltényi.com

**South Korea**

Miltényi Biotec Korea Co., Ltd.  
Arigi Bldg, 8F  
562 Nonhyeon-ro  
Gangnam-gu  
Seoul 06136, South Korea  
Phone +82 2 555 1988  
Fax +82 2 555 8890  
macskr@miltényi.com

**Spain**

Miltényi Biotec S.L.  
C/Luis Buñuel 2  
Ciudad de la Imagen  
28223 Pozuelo de Alarcón (Madrid)  
Spain  
Phone +34 91 512 12 90  
Fax +34 91 512 12 91  
macses@miltényi.com

**Switzerland**

Miltényi Biotec Swiss AG  
Gibelinstrasse 27  
4500 Solothurn  
Switzerland  
Phone +41 32 623 08 47  
Fax +49 2204 85197  
macsch@miltényi.com

**United Kingdom**

Miltényi Biotec Ltd.  
Almac House, Church Lane  
Bisley, Surrey GU24 9DR, UK  
Phone +44 1483 799 800  
Fax +44 1483 799 811  
macsuk@miltényi.com

**www.miltényibiotec.com**

Miltényi Biotec provides products and services worldwide. Visit [www.miltényibiotec.com/local](http://www.miltényibiotec.com/local) to find your nearest Miltényi Biotec contact.

Unless otherwise specifically indicated, Miltényi Biotec products and services are for research use only and not for therapeutic or diagnostic use.

The CliniMACS System components, including Reagents, Tubing Sets, Instruments, and PBS/EDTA Buffer, are designed, manufactured and tested under a quality system certified to ISO 13485.

In the EU, the CliniMACS System components are available as CE-marked medical devices for their respective intended use, unless otherwise stated. The CliniMACS Reagents and Biotin Conjugates are intended for *in vitro* use only and are not designated for therapeutic use or direct infusion into patients. The CliniMACS Reagents in combination with the CliniMACS System are intended to separate human cells. Miltényi Biotec as the manufacturer of the CliniMACS System does not give any recommendations regarding the use of separated cells for therapeutic purposes and does not make any claims regarding a clinical benefit. For the manufacturing and use of target cells in humans, the national legislation and regulations – e.g. for the EU the Directive 2004/23/EC (“human tissues and cells”), or the Directive 2002/98/EC (“human blood and blood components”) – must be followed. Thus, any clinical application of the target cells is exclusively within the responsibility of the user of a CliniMACS System.

In the US, the CliniMACS CD34 Reagent System, including the CliniMACS Plus Instrument, CliniMACS CD34 Reagent, CliniMACS Tubing Sets TS and LS, and the CliniMACS PBS/EDTA Buffer, is FDA approved as a Humanitarian Use Device (HUD), authorized by U.S. Federal law for use in the treatment of patients with acute myeloid leukemia (AML) in first complete remission. The effectiveness of the device for this indication has not been demonstrated. Other products of the CliniMACS Product Line are available for use only under an approved Investigational New Drug (IND) application, Investigational Device Exemption (IDE) or FDA approval.

CliniMACS GMP MicroBeads are for research use and *ex vivo* cell processing only.

CliniMACS MicroBeads are for research use only and not for human therapeutic or diagnostic use.

MACS GMP Products are for research use and *ex vivo* cell culture processing only, and are not intended for human *in vivo* applications. For regulatory status in the USA, please contact your local representative. MACS GMP Products are manufactured and tested under a quality system certified to ISO 13485 and are in compliance with relevant GMP guidelines. They are designed following the recommendations of USP <1043> on ancillary materials.

CentriCult, CliniMACS, CliniMACS Prodigy, MACS, PepTivator, TexMACS, and the Miltényi Biotec logo are registered trademarks or trademarks of Miltényi Biotec and/or its affiliates in various countries worldwide. Copyright © 2022 Miltényi Biotec and/or its affiliates. All rights reserved.